

PAWR Antibody



PACO11080

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAWR

Uniprot

Q96IZ0

Synonyms:

PRKC, apoptosis, WT1, regulator; PAWR; PAR4; Par-4 ;

Immunogen:

Human PAWR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A