

# PAX6 Antibody



PACO11083

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PAX6

**Uniprot**

P26367

**Synonyms:**

paired box 6; PAX6; AN; AN2; D11S812E; MGC17209; MGDA; WAGR ;

**Immunogen:**

Human PAX6.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A