

PCBP2 Antibody



PACO11090

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PCBP2

Uniprot

Q15366

Synonyms:

poly(rC) binding protein 2; PCBP2; HNRPE2; MGC110998; hnRNP-E2 ;

Immunogen:

Human PCBP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A