

# PCK1 Antibody



PACO11105

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCK1

**Uniprot**

P35558

**Synonyms:**

phosphoenolpyruvate carboxykinase 1 (soluble); PCK1; MGC22652; PEPCK-C; PEPCK1; PEPCKC ;

**Immunogen:**

Human PCK1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A