

PCOLCE2 Antibody



PACO11112

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCOLCE2

Uniprot

Q9UKZ9

Synonyms:

procollagen C-endopeptidase enhancer 2; PCOLCE2; PCPE2 ;

Immunogen:

Human PCOLCE2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A