

PCSK5 Antibody



PACO11114

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCSK5

Uniprot

Q92824

Synonyms:

proprotein convertase subtilisin/kexin type 5; PCSK5; PC5; PC6; PC6A; SPC6 ;

Immunogen:

Human PCSK5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A