

PCSK7 Antibody



PACO11115

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCSK7

Uniprot

Q16549

Synonyms:

proprotein convertase subtilisin/kexin type 7; PCSK7; LPC; PC7; PC8; SPC7 ;

Immunogen:

Human PCSK7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A