

CDK18 Antibody



PACO11118

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CDK18

Uniprot

Q07002

Synonyms:

PCTAIRE protein kinase 3; PCK3; PCTAIRE; PCTAIRE3 ;

Immunogen:

Human PCK3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A