

PCTP Antibody



PACO11119

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PCTP

Uniprot

Q9UKL6

Synonyms:

phosphatidylcholine transfer protein; PCTP; STARD2 ;

Immunogen:

Human PCTP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A