

PCYT1B Antibody



PACO11121

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PCYT1B

Uniprot

Q9Y5K3

Synonyms:

phosphate cytidylyltransferase 1, choline, beta; PCYT1B; CCT-beta; CTB ;

Immunogen:

Human PCYT1B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A