

# PDCD10 Antibody



PACO11125

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PDCD10

**Uniprot**

Q9BUL8

**Synonyms:**

programmed cell death 10; PDCD10; CCM3; MGC1212; MGC24477; TFAR15 ;

**Immunogen:**

Human PDCD10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A