

PDCD6 Antibody



PACO11129

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PDCD6

Uniprot

O75340

Synonyms:

programmed cell death 6; PDCD6; ALG-2; FLJ46208; MGC111017; MGC119050; MGC9123; PEF1B ;

Immunogen:

Human PDCD6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A