

PDK4 Antibody



PACO11145

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PDK4

Uniprot

Q16654

Synonyms:

pyruvate dehydrogenase kinase, isozyme 4; PDK4; FLJ40832 ;

Immunogen:

Human PDK4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A