

# PDS5A Antibody



PACO11151

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PDS5A

**Uniprot**

Q29RF7

**Synonyms:**

PDS5, regulator of cohesion maintenance, homolog A (S. cerevisiae); PDS5A; DKFZp686B19246; FLJ41012; KIAA0648; MGC131948; MGC161503; PIG54; SCC-112 ;

**Immunogen:**

Human PDS5A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A