

# PDS5B Antibody



PACO11152

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PDS5B

**Uniprot**

Q9NTI5

**Synonyms:**

PDS5, regulator of cohesion maintenance, homolog B (S. cerevisiae); PDS5B; APRIN; AS3; CG008; FLJ23236; KIAA0979; RP1-267P19.1 ;

**Immunogen:**

Human PDS5B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A