

PDSS2 Antibody



PACO11153

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PDSS2

Uniprot

Q86YH6

Synonyms:

prenyl (decaprenyl) diphosphate synthase, subunit 2; PDSS2; C6orf210; DLP1; bA59I9.3; hDLP1 ;

Immunogen:

Human PDSS2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A