

PACO11155

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PDXK

Uniprot

O00764

Synonyms:

pyridoxal (pyridoxine, vitamin B6) kinase; PDXK; C21orf124; C21orf97; DKFZp566A071; FLJ31940; FLJ37311; MGC15873; MGC31754; MGC52346; PKH; PNK; PRED79 ;

Immunogen:

Human PDXK.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A