

# PDZD2 Antibody



PACO11156

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PDZD2

**Uniprot**

O15018

**Synonyms:**

PDZ domain containing 2; PDZD2; AIPC; KIAA0300; PAPIN; PDZK3; PIN1 ;

**Immunogen:**

Human PDZD2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A