

PDZD3 Antibody



PACO11157

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PDZD3

Uniprot

Q86UT5

Synonyms:

PDZ domain containing 3; PDZD3; FLJ22756; IKEPP; PDZK2 ;

Immunogen:

Human PDZD3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A