

PEBP1 Antibody



PACO11159

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PEBP1

Uniprot

P30086

Synonyms:

phosphatidylethanolamine binding protein 1; PEBP1; HCNP; PBP; PEBP; RKIP ;

Immunogen:

Human PEBP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A