

PEF1 Antibody



PACO11163

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PEF1

Uniprot

Q9UBV8

Synonyms:

penta-EF-hand domain containing 1; PEF1; PEF1A; PEFLIN ;

Immunogen:

Human PEF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A