

PEG10 Antibody



PACO11164

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PEG10

Uniprot

Q86TG7

Synonyms:

paternally expressed 10; PEG10; Edr; HB-1; KIAA1051; MEF3L; Mar2; Mart2; RGAG3 ;

Immunogen:

Human PEG10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A