

PELO Antibody



PACO11166

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PELO

Uniprot

Q9BRX2

Synonyms:

pelota homolog (Drosophila); PELO; CGI-17; PRO1770 ;

Immunogen:

Human PELO.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A