

PEPD Antibody



PACO11167

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PEPD

Uniprot

P12955

Synonyms:

peptidase D; PEPD; MGC10905; PROLIDASE ;

Immunogen:

Human PEPD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A