

PER1 Antibody



PACO11168

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

PER1

Uniprot

O15534

Synonyms:

period homolog 1 (Drosophila); PER1; MGC88021; PER; RIGUI; hPER ;

Immunogen:

Human PER1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A