

PER2 Antibody



PACO11169

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PER2

Uniprot

O15055

Synonyms:

period homolog 2 (Drosophila); PER2; FASPS; KIAA0347 ;

Immunogen:

Human PER2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A