

PEX1 Antibody



PACO11172

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PEX1

Uniprot

O43933

Synonyms:

peroxisomal biogenesis factor 1; PEX1; ZWS1 ;

Immunogen:

Human PEX1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A