

PEX11G Antibody



PACO11173

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PEX11G

Uniprot

Q96HA9

Synonyms:

peroxisomal biogenesis factor 11 gamma; PEX11G; MGC4281 ;

Immunogen:

Human PEX11G.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A