

PFKM Antibody



PACO11186

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PFKM

Uniprot

P08237

Synonyms:

phosphofructokinase, muscle; PFKM; GSD7; MGC8699; PFK-1; PFK-M; PFKX ;

Immunogen:

Human PFKM.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A