

PGM2 Antibody



PACO11200

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PGM2

Uniprot

Q96G03

Synonyms:

phosphoglucomutase 2; PGM2; FLJ10983; MSTP006 ;

Immunogen:

Human PGM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A