

PGM2L1 Antibody



PACO11201

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PGM2L1

Uniprot

Q6PCE3

Synonyms:

phosphoglucomutase 2-like 1; PGM2L1; BM32A; FLJ32029 ;

Immunogen:

Human PGM2L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A