

PGPEP1 Antibody



PACO11202

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PGPEP1

Uniprot

Q9NXJ5

Synonyms:

pyroglutamyl-peptidase I; PGPEP1; FLJ20208; MGC10812; PGP; PGP-I; PGPI; Pcp ;

Immunogen:

Human PGPEP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A