

PGS1 Antibody



PACO11204

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PGS1

Uniprot

Q32NB8

Synonyms:

phosphatidylglycerophosphate synthase 1; PGS1; DKFZp762M186; MGC131960 ;

Immunogen:

Human PGS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A