

PHF11 Antibody



PACO11209

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PHF11

Uniprot

Q9UIL8

Synonyms:

PHD finger protein 11; PHF11; APY; BCAP; IGEL; IGER; IGER; NY-REN-34; NYREN34 ;

Immunogen:

Human PHF11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A