

PHF6 Antibody



PACO11213

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PHF6

Uniprot

Q8IWS0

Synonyms:

PHD finger protein 6; PHF6; BORJ; MGC14797 ;

Immunogen:

Human PHF6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A