

PHKB Antibody



PACO11216

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PHKB

Uniprot

Q93100

Synonyms:

phosphorylase kinase, beta; PHKB; DKFZp781E15103; FLJ41698 ;

Immunogen:

Human PHKB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A