

PHKG1 Antibody



PACO11217

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PHKG1

Uniprot

Q16816

Synonyms:

phosphorylase kinase, gamma 1 (muscle); PHKG1; PHKG ;

Immunogen:

Human PHKG1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A