

PHKG2 Antibody



PACO11218

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PHKG2

Uniprot

P15735

Synonyms:

phosphorylase kinase, gamma 2 (testis); PHKG2; GSD9C ;

Immunogen:

Human PHKG2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A