

PHLDA2 Antibody



PACO11220

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PHLDA2

Uniprot

Q53GA4

Synonyms:

pleckstrin homology-like domain, family A, member 2; PHLDA2; BRW1C; BWR1C; HLDA2; IPL; TSSC3 ;

Immunogen:

Human PHLDA2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A