

PHYHD1 Antibody



PACO11223

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PHYHD1

Uniprot

Q5SRE7

Synonyms:

phytanoyl-CoA dioxygenase domain containing 1; PHYHD1; MGC16638 ;

Immunogen:

Human PHYHD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A