

# PHYHIP Antibody



PACO11224

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PHYHIP

**Uniprot**

Q92561

**Synonyms:**

phytanoyl-CoA 2-hydroxylase interacting protein; PHYHIP; DYRK1AP3; KIAA0273; PAHX-AP; PAHXAP1 ;

**Immunogen:**

Human PHYHIP.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A