

PI16 Antibody



PACO11226

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PI16

Uniprot

Q6UXB8

Synonyms:

peptidase inhibitor 16; PI16; CRISP9; DKFZp586B1817; MGC45378; MSMBBP; PSPBP ;

Immunogen:

Human PI16.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A