

PI4K2A Antibody



PACO11227

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PI4K2A

Uniprot

Q9BTU6

Synonyms:

phosphatidylinositol 4-kinase type 2 alpha; PI4K2A; DKFZp761G1923; PI4KII; PIK42A; RP11-548K23.6 ;

Immunogen:

Human PI4K2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A