

PI4KB Antibody



PACO11230

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PI4KB

Uniprot

Q9UBF8

Synonyms:

phosphatidylinositol 4-kinase, catalytic, beta; PI4KB; PI4K-BETA; PI4KIIIbeta; PI4Kbeta; PIK4CB; pi4K92 ;

Immunogen:

Human PI4KB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A