

# PIBF1 Antibody



PACO11234

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PIBF1

**Uniprot**

Q8WXW3

**Synonyms:**

progesterone immunomodulatory binding factor 1; PIBF1; C13orf24; KIAA1008; PIBF; RP11-505F3.1 ;

**Immunogen:**

Human PIBF1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A