

PIK3IP1 Antibody



PACO11247

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PIK3IP1

Uniprot

Q96FE7

Synonyms:

phosphoinositide-3-kinase interacting protein 1; PIK3IP1; HGFL; MGC17330; hHGFL(S) ;

Immunogen:

Human PIK3IP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A