

PIK3R3 Antibody



PACO11249

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PIK3R3

Uniprot

Q92569

Synonyms:

phosphoinositide-3-kinase, regulatory subunit 3 (gamma); PIK3R3; DKFZp686P05226; FLJ41892; p55; p55-gamma ;

Immunogen:

Human PIK3R3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A