

PIK3R5 Antibody



PACO11251

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PIK3R5

Uniprot

Q8WYR1

Synonyms:

phosphoinositide-3-kinase, regulatory subunit 5; PIK3R5; F730038I15Rik; FOAP-2; P101-PI3K; p101 ;

Immunogen:

Human PIK3R5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A