

PIN1 Antibody



PACO11252

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PIN1

Uniprot

Q13526

Synonyms:

peptidylprolyl cis/trans isomerase, NIMA-interacting 1; PIN1; DOD; UBL5 ;

Immunogen:

Human PIN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A