

PITPNB Antibody



PACO11262

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PITPNB

Uniprot

P48739

Synonyms:

phosphatidylinositol transfer protein, beta; PITPNB; PI-TP-beta; PtdInsTP; VIB1B ;

Immunogen:

Human PITPNB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A